

PBUI09/02.08.01/sx

15

Abstract

A sample vessel with a sample to be evaporated is held in a holding device. With the holding device the sample vessel may be placed into an evaporator.

The receiving opening (3) of the holding device (1) comprises a cooling section (4) which is provided with cooling means. Thanks to the cooling means an exactly definable residual volume (R) in the sample vessel (2) is not evaporated.

02-08-01 15:09 VON -Hepp Wenger & Ryffel AG, Wil +41-71-9139556 T-330 P.018/055 F-557