

02-08-01 15:08 VON -Hepp Wenger &amp; Ryffel AG, Wil

+41-71-9139556

T-330 P.018/055 F-557

PBUL09/02.08.01/ax

15

### Abstract

A sample vessel with a sample to be evaporated is held in a holding device. With the holding device the sample vessel may be placed into an evaporator.

The receiving opening (3) of the holding device (1) comprises a cooling section (4) which is provided with cooling means. Thanks to the cooling means an exactly definable residual volume (R) in the sample vessel (2) is not evaporated.

02-08-01 15:08 VON -Hepp Wenger & Ryffel AG, Wil  
+41-71-9139556  
T-330 P.018/055 F-557